



ULT-2020 and Probe Adapter Validated by Philips Ultrasound

After many months of conducting an exhaustive evaluation of the BC Group's Ultrasound Leakage Tester and Probe Adapter, Philips Ultrasound has issued the following statement...

[Read more...](#)



April 2016
www.bcgrouppintl.com

MD EXPO EXHIBIT HALL HOURS

Friday, April 22:
1:00 PM - 4:00 PM
Saturday, April 23:
12:30 PM - 3:30 PM

Booth 502
(Near the Entrance)

[Visit the BC Marketplace](#)

Download Your Free Pass to MD Expo 2016

On April 21-23, 2016, the MD Expo will be held at the Fairmont Hotel in Dallas, Texas. You can download a free pass to the show, which includes education courses approved for CE Credit, six hours of exhibit hall time with over a hundred vendors, multiple opportunities for networking and more!

[Read more...](#)

SA-2500 Safety Analyzer

The new SA-2500 is a portable, versatile and accurate safety analyzer from BC Group. The SA-2500 is unique in that it offers a Remote Mode of operation. Users can control the safety analyzer directly from their PC. All measurements can be integrated into the individual operator interface at the PC and to define test sequences.



SA-2500

[Read more...](#)



NIBP-1040 Kit

For complete patient simulation, BC Group offers the NIBP-1040 Multi-Parameter Patient Simulator. The NIBP-1040 handles a wide variety of functions, including non-invasive blood pressure, invasive blood pressure, ECG, arrhythmias, respiration, leak detection, fetal/maternal and cardiac output. Each feature can be accessed through the graphic display. Users can view numerical data or graphical representation of data.



NIBP-1040 Kit

[Read more...](#)



TSI Certifier FA Plus

2015-2016 Catalogs



Ventilation Testing with the PFC-3000 Benchtop Flow Analyzer



PFC-3000

The PFC-3000 series of ventilator testers measures flow, pressure, temperature, humidity and O2 concentrations bidirectionally. BC Group is the only factory authorized service center for the PFC-3000 series and IMT products in North America.

[Read more...](#)

Contact Us...

Phone: 1 (314) 638-3800
Hours: 7:00AM - 5:30PM
Email: sales@bcgrouppintl.com

BC Group International, Inc.
3081 Elm Point Industrial Dr | St. Charles, MO 63301-4333 USA
sales@bcgrouppintl.com | service@bcgrouppintl.com | 1-888-223-6763 | 1-314-638-3800



[Subscribe](#) to our email list