



Ultrasound Leakage Testing

The Model ULT-2000 Series is a family of microprocessor-based Ultrasound Transducer Leakage Testers. The ULT-2010 measures both the conductivity of the cleaning medium and the leakage current of the ultrasound transducer. The ULT-2020 offers the same features of the ULT-2010, plus a Meter Mode and data logging.



[Read More...](#)

Download the ForceTriad™ Verison 4.0 Autosequence Upgrade for the ESU-2400H



The newest ESU-2400H update includes bug fixes and compatibility with the ForceTriad up to version 4.0. In order to perform autosequences on the ForceTriad™, users will need the [ForceTriad PM Kit](#) (BC20-00131).

[Read More...](#)

July 2017

www.bcgrouptl.com

Download the ULT Comparison Chart



[View BC Group's Upcoming Events](#)

Contact Us...

Phone: 1 (314) 638-3800

Hours: 8:00AM - 5:00PM CST

Email: sales@bcgrouptl.com

BC Group International Inc.
3081 Elm Point Industrial Dr. | St. Charles, MO 63301-4333 USA
sales@bcgrouptl.com | service@bcgrouptl.com | 1-888-223-6763 | 1-314-638-3800



[Subscribe](#) to our email list