



Riken FI-8000P Portable Gas Indicator



When needing to test anesthetic vaporizer output concentration, the Riken FI-8000P is the most capable gas indicator available. It is the latest development in the line of Riken gas indicators based on the optical interferometer principle. The Riken FI-8000P makes all calculations and adjustments for range, anesthetic and conversion factors automatically.

[Read More...](#)

Ultrasound Leakage Testing

The Model ULT-2000 Series is a family of microprocessor-based Ultrasound Transducer Leakage Testers. The ULT-2010 measures both the conductivity of the cleaning medium and the leakage current of the ultrasound transducer. The ULT-2020 offers the same features of the ULT-2010, plus a Meter Mode and data logging.



[Read More...](#)



August 2017

www.bcgrouptl.com

Shop NEW RTI Imaging



Test Equipment

[Download Your FREE PASS to MD Expo](#)

MD EXPO

Contact Us...

Phone: 1 (314) 638-3800

Hours: 8:00AM - 5:00PM CST

Email: sales@bcgrouptl.com

BC Group International Inc.
3081 Elm Point Industrial Dr. | St. Charles, MO 63301-4333 USA
sales@bcgrouptl.com | service@bcgrouptl.com | 1-888-223-6763 | 1-314-638-3800



[Subscribe](#) to our email list