



Think Outside the Box with the IPA-3400



The new IPA-3400 Infusion Pump Analyzer is designed to minimize downtime. It is the only multiple channel device to provide OUTSIDE THE BOX repair and calibration by allowing the user to easily remove the modules separately. You can still perform testing with the remaining modules providing a continuous uptime operation.

[Read More...](#)

Now Shipping: NEW Anesthetic Agent Analyzer

The AA-8000 is the latest in portable anesthetic agent analyzers. It's designed to meet the demand for a small, easy to use unit with high reliability and accuracy. Utilizing proprietary state-of-the-art digital NDIR (NonDispersive Infrared) Technology, the AA-8000 provides a low cost, high function, microprocessorbased analyzer that is simple to operate while maintaining high performance and accuracy.



The AA-8000 also comes with a 12 point NIST Traceable Calibration Certificate with Data! [View a sample](#) of the calibraton data to see what is covered

[Read More...](#)

November 2018

www.bcgrouptl.com

Visit BC Group at



Booth 107

**Pearson Current
Transformer Calibration
Now Available from BC
Group**



Contact Us...

Phone: 1 (314) 638-3800

Hours: 8:00AM - 5:00PM CST

Email: sales@bcgrouptl.com