



IPA-3400 Infusion Pump Analyzer Receives US Patent

The IPA-3400 Infusion Pump Analyzer has received US Patent No. 10,100,828. The IPA-3400 features a new, cutting-edge, dual syringe stepper motor driven system that provides a more realistic flow path for the infusion device under test. Also, independent stepper motor control of the custom designed ceramic valving allows the system to run not only more quietly and more smoothly, but it also allows for a bidirectional powered fluid flow for use in the built-in cleaning cycle.



[Learn More...](#)

Think Outside the Box with the IPA-3400



The new IPA-3400 Infusion Pump Analyzer is designed to minimize downtime. It is the only multiple channel device to provide OUTSIDE THE BOX repair and calibration by allowing the user to easily remove the modules separately. You can still perform testing with the remaining modules providing a continuous uptime operation.

[Read More...](#)

December 2018

www.bcgrouptl.com

For Hands On ESU Trainging Visit



BIOMED
SUPPORT SERVICES
Simplifying the Complex

Shop Atrix Vacuums



Contact Us...

Phone: 1 (314) 638-3800

Hours: 8:00AM - 5:00PM CST

Email: sales@bcgrouptl.com

BC Group International Inc.
3081 Elm Point Industrial Dr. | St. Charles, MO 63301-4333 USA
sales@bcgrouptl.com | service@bcgrouptl.com | 1-888-223-6763 | 1-314-638-3800

emma[®]