




## Biomed Support Services Is Bringing ESU Training to Your Area

BC Group and Biomed Support Services are pleased to announce a series of comprehensive, hands-on training courses for 2019. These courses are available to anyone interested in obtaining the very best information on a variety of Valleylab electrosurgical generators. The instructor, Dale Munson, is back with his simple-as-possible training materials and a passion for helping biomed.

[Read More...](#)

## IPA-3400 Infusion Pump Analyzer Receives US Patent

The IPA-3400 Infusion Pump Analyzer has received US Patent No. 10,100,828. The IPA-3400 features a new, cutting-edge, dual syringe stepper motor driven system that provides a more realistic flow path for the infusion device under test. Also, independent stepper motor control of the custom designed ceramic valving allows the system to run not only more quietly and more smoothly, but it also allows for a bidirectional powered fluid flow for use in the built-in cleaning cycle.



[Learn More...](#)


January 2019

[www.bcgrouptl.com](http://www.bcgrouptl.com)

Download your **FREE PASS**

# MD EXPO

Shop **ESU Analyzers**



**Contact Us...**

**Phone:** 1 (314) 638-3800

**Hours:** 8:00AM - 5:00PM CST

**Email:** [sales@bcgrouptl.com](mailto:sales@bcgrouptl.com)

BC Group International Inc.  
3081 Elm Point Industrial Dr. | St. Charles, MO 63301-4333 USA  
[sales@bcgrouptl.com](mailto:sales@bcgrouptl.com) | [service@bcgrouptl.com](mailto:service@bcgrouptl.com) | 1-888-223-6763 | 1-314-638-3800

emma