



## Happy Holidays from BC Group



Happy Holidays from all of us at BC Group! Due to the upcoming holidays, our office will be closed **December 24th-25th** and **January 1st**. We hope you stay safe and enjoy the holiday season.

## Shop for any Last Minute Needs



With the season of giving nearly here be sure to visit our store for any needs you might have! We offer many products from ventilator analyzers to patient simulators and infusion pump analyzers. BC Group also offers a wide variety of other products and brands all at the click of a button and they can all be found on our webstore!

BC Group International Inc.  
3081 Elm Point Industrial Dr. | St. Charles, MO 63301-4333 USA  
[sales@bcgroupintl.com](mailto:sales@bcgroupintl.com) | [service@bcgroupintl.com](mailto:service@bcgroupintl.com) | 1-888-223-6763 | 1-314-638-3800

December 2020  
[www.bcgrouptnl.com](http://www.bcgrouptnl.com)

## Get the NEW FSX-1101



The FSX-1101 is a digital SPO2 simulator with many features, including a **LIFETIME Warranty\***. It is small and light weight, yet very powerful. It has a bright, 2.4" color touchscreen that, coupled with an intuitive menu structure, allows easy access and viewing of all setup and function parameters.

### Features - FSX-1101:

- Small, Handheld, Lightweight, Rugged
- Intuitive User Interface
- Color Touch Screen
- Preset & Manual Mode Operation
- 6 Saturation Levels from 80-99%
- Heart Rates from 30-240 BPM
- Adjustable Perfusion Index
- Pulse Synchronized with NIBP-1000, PS-2100 or PS-2200 Series

\* LIFETIME Limited Warranty requires product registration and return to Factory for annual calibration.

### Contact Us...

**Phone:** 1 (314) 638-3800  
**Hours:** 8:00AM - 5:00PM CST  
**Email:** [sales@bcgroupintl.com](mailto:sales@bcgroupintl.com)

emma