



Teach and Improve Patient Positioning with CIRS Products Sold Through BC Group



An ideal addition to any imaging department or Radiologic Technology training program, the CIRS products can assist in the monitoring, training and improvement of parameters and protocols common to most pediatric imaging procedures. The phantom is designed to assist technical and clinical staff in the selection, monitoring, training and verification of scanning parameters common to most radiological procedures requiring fine anatomical details.

Benefits:

- Facilitates hands on training and feedback in an exciting memorable way
- Demonstrates the need for specialized views with certain pathologies
- Allows for effective demonstration of changes in techniques
- Promotes sensitivity and proficiency in patient handling and positioning
- Effectively evaluate trainee performance

Keep Up To Date on Changes to Ultrasound Leakage Testing for Your TEE Probes



Did you know that if you are using some GE TEE probes or Sonosite TEE probes that you need to test the leakage of the probe and for possible leakage at the control knob? To learn more about how to test the probes and any special requirements for some Manufacture's, go to the following link and view the procedures and videos on how to perform these simples tests.

[Learn More...](#)

February 2020
www.bcgrouptl.com

Shop Anesthetic Agent Analyzer-8000



The AA-8000 is the latest in portable anesthetic agent analyzers. Utilizing proprietary state-of-the-art digital NDIR (Non-Dispersive Infrared) Technology it provides a low cost, high function, microprocessor-based analyzer that is simple to operate while maintaining high performance and accuracy.

Shop Non-Invasive Blood Pressure & Multiparamator-1040 Kit



Test 13 parameters with one portable device. Most cost effective multiparameter simulator on the market.

Contact Us...

Phone: 1 (314) 638-3800

Hours: 8:00AM - 5:00PM CST

Email: sales@bcgrouptl.com

BC Group International Inc.
3081 Elm Point Industrial Dr. | St. Charles, MO 63301-4333 USA
sales@bcgrouptl.com | service@bcgrouptl.com | 1-888-223-6763 | 1-314-638-3800

emma®