




## COVID-19 Calibration & Repair Safety Precautions

At BC Group International, Inc., there is no higher priority than the health and safety of our employees and community. There is no time more important than now in delivering lifesaving, quality products and services to those in need.



Our Service team continues to provide Critical Essential Services for our customers. We are taking important steps to maintain our services while safeguarding your and our employees' health and safety.

[Learn More...](#)

## Check out the DA-2006 Series



The DA-2006 Series Defibrillator Analyzers provides Variable Loads used when testing Defibrillators to assure the proper electrical current is delivered to the heart, per IEC 60601-2-4 and AAMI DF80 standards with automatic detection of loads. All models have a built in 50  $\Omega$  human body simulation load as well as 12 lead ECG with arrhythmias and performance waveforms. Additionally, they have a Parallel Printer port, a RS-232 serial communications port, oscilloscope output, high-level ECG output, as well as provision for an AC adaptor.

July 2020

[www.bcgrouptl.com](http://www.bcgrouptl.com)

## Shop Ventilator Testers



The PFC-3000 series of ventilator testers measures flow, pressure, temperature, humidity and O<sub>2</sub> concentrations bidirectionally. It combines a simple, intuitive multilingual user interface with the highest precision tester.

## Test Your Infusion Pump Easily



The IPA-3400 is a high accuracy, easy to use system that incorporates full touch screen control of all processes without the use of old fashioned buttons and knobs.

### Contact Us...

Phone: 1 (314) 638-3800

Hours: 8:00AM - 5:00PM CST

Email: [sales@bcgrouptl.com](mailto:sales@bcgrouptl.com)

BC Group International Inc.  
3081 Elm Point Industrial Dr. | St. Charles, MO 63301-4333 USA  
[sales@bcgrouptl.com](mailto:sales@bcgrouptl.com) | [service@bcgrouptl.com](mailto:service@bcgrouptl.com) | 1-888-223-6763 | 1-314-638-3800

emma