



Get Your New Ventilator Analyzer Fast

In Stock and Ready to Ship



BC Group now has several PFC-3000A's and PFC3000L's in stock, which means any order placed will be shipped to you quickly. The PFC-3000 series of ventilator testers measures flow, pressure, temperature, humidity and O2 concentrations bidirectionally. The one-of-a-kind adult, pediatric and high frequency ventilation measuring modes make the PFC-3000 the ideal calibration tool for all ventilators, anesthesia machines and spirometers. The PFC-3000 distinguishes itself from other calibration tools by combining a simple, intuitive multilingual user interface with the highest precision tester. Don't miss out on the opportunity to get your vital test equipment in just a few weeks.

COVID-19 Calibration & Repair Safety Precautions

At BC Group International, Inc., there is no higher priority than the health and safety of our employees and community. There is no time more important than now in delivering lifesaving, quality products and services to those in need.



Our Service team continues to provide Critical Essential Services for our customers. We are taking important steps to maintain our services while safeguarding your and our employees' health and safety.

[Learn More...](#)

August 2020
www.bcgrouptl.com

[Shop NIBP Series](#)



BC Biomedical NIBP Simulators offer the best bang for your buck in today's market. The NIBP Series is a Microprocessor based, High Precision Non-Invasive Blood Pressure (NIBP) Simulators.

[Explore Anesthetic Agent Analyzer-8000](#)



The AA-8000 is the latest in portable anesthetic agent analyzers. Utilizing proprietary state-of-the-art digital NDIR (Non-Dispersive Infrared) Technology it provides a low cost, high function, microprocessor-based analyzer that is simple to operate while maintaining high performance and accuracy.

Contact Us...

Phone: 1 (314) 638-3800
 Hours: 8:00AM - 5:00PM CST
 Email: sales@bcgroupintl.com

BC Group International Inc.
 3081 Elm Point Industrial Dr. | St. Charles, MO 63301-4333 USA
sales@bcgroupintl.com | service@bcgroupintl.com | 1-888-223-6763 | 1-314-638-3800

emma®